

## **PVT500 the self-sufficient test system for pvT-measurement – A system comparison**

pvT-measurements are used to determine the thermodynamic behaviour of plastics, elastomers, thermosets and liquids. Here, the specific volume is determined as a function of the variables of pressure and temperature. pvT-measurements are essential for the following applications:

- Molded part design in injection molding
- Material selection at existing injection molds
- Shrinking processes in extrusion
- Shrinking processes when pressing moldings
- Determination of shrinkage caused by vulcanization or linkage
- Design of expansion tanks for liquids

**PVT500** operates according to **ISO17744**, which distinguishes between isothermal and isobaric measuring method. During isothermal measuring, the sample is exposed to different pressures at a constant temperature. The isobaric measurement method is carried out at constant pressure, wherein the sample is cooled at a defined cooling rate.

The isothermal measuring method is still one of the most common applications, however the isobaric measuring method describes the process significantly better. This is especially the case in injection molding because the molding cools down at constant pressures (upto sealing- or freezing point of the sprue system). Such a handling process is usually used at current tool designs.

### **Features of PVT500 in comparison to Add-on option\* „Capillary Rheometer“**

The new **PVT500** offers several advantages in comparison to the optional available Add-on for Capillary Rheometer:

- Test barrel with just 9,5mm diameter allows a faster cooling of the sample with better temperature homogeneity.
- Increase of cooling rate to **40K/min** – compared to previous possible 25K/min at the Add-on option.
- Maximum pressure up to 2500 bar
- Less amount of sample necessary
- Stand-alone device for higher test capacity
- Higher accuracy compared to Add-on option
- Smaller sample provides more homogeneous stress state in solid state range, especially with amorphous plastics

\* The **Add-on option "PVT"**, is one of various extensions which is used successfully with GOETTFERT Capillary Rheometer.

## Comparison between PVT500 and Add-on option

Test: Isobaric – cooling rate 2,5K/min  
 Test material: Polycarbonate  
 RHEOGRAPH: Liquid cooling via oil  
 PVT500: Cooling via compressed air

The comparative tests were done with **PVT500** and **RHEOGRAPH 75** with integrated **Add-on option PVT** in isobaric measuring mode at a cooling rate of 2,5 K/min.




Figure 1: PVT500




Figure 2: Add-on (RHEOGRAPH 75)

Both diagrams show the same course. Though the Isobars have a more linear course with PVT 500, so this data can be approximated better with Tait-Model. The data will be approximated with the Tait-model and the 1bar pressure line will be extrapolated for evaluation (see following explanation). There are two measuring points of this 1bar-Isobar (orange arrows figure 1 and 2) which can easily be compared with other methods:

Density at room temperature (Determination by **buoyancy method** 23°C) and density at melt temperature (Determination by Melt Flow Indexer 300°C).

**Table 1: Comparison of density values for polycarbonate at RT and 300°C**

Test temperature	Density determination	PVT500 (Tait)	RHEOGRAPH75 (Tait)
23°C	1,2 g/cm <sup>3</sup>	1,189g/cm <sup>3</sup> (-1%)	1,222 g/cm <sup>3</sup> (1,8%)
300°C	1,065g/cm <sup>3</sup>	1,077g/cm <sup>3</sup> (1,1%)	1,104 g/cm <sup>3</sup> (3,6%)

The comparison of the data indicates, that the approximated data of **PVT500** shows a 1% deviation from the density of other methods. On the other hand, the approximated data with the Add-on of Capillary Rheometer shows an almost three times higher deviation.

Similar tests were performed with different polypropylene.

**Table 2: Comparison of density values for PP at RT/30°C and 230°C**

Test temperature	Density determination	PVT500 (Tait)	RHEOGRAPH75 (Tait)
23°C	0,903 g/cm <sup>3</sup>	0,900 g/cm <sup>3</sup> (-0,4%)	0,909 g/cm <sup>3</sup> (0,7%)
230°C	0,74 g/cm <sup>3</sup>	0,741 g/cm <sup>3</sup> (0,1%)	0,746 g/cm <sup>3</sup> (0,8%)

Here you can see the same effect. PVT500 exceeds the already high accuracy of the Add-on option (**RHEOGRAPH 75**).



### Example data of different plastics

The following diagrams contain some sample data for PEEK and PEI. The data was determined with a cooling rate of 5K/min in isobaric measuring mode.

The data of PEEK show the typical behaviour for a semi-crystalline material with distinct transition zone at the melting point. The data of PEI show the typical amorphous behaviour with smooth transition in melting range.

The measurement of PEI with amorphous behaviour shows, with its better linear behaviour of the isobars, an improvement in comparison to Add-on option Capillary Rheometer.

The measured data were additionally approximated with the Tait model. On this basis, the Isobar at ambient pressure (p=0bar) was determined. The data is accurately described and summarized by the model.



## Modeling the measurement data via Tait-model

The main application for modelling of pVT measurement data via the Tait-model is in simulation of injection molded components.

Because of the modelling, one is no longer limited to discretely measured Isobars. One can also rely on interpolated model data depending on which operating point is set on the injection molding machine and how the process continues in the tool.

Such an operation is shown as example in the following diagram.

After injection phase, the holding-pressure phase takes place up until the sealing point. Then a further cooling process follows up until the tool opens.





Figure 3: Tait-Modell

The sealing point is very important in this process. It can be determined directly from the Tait-model. The model describes the sectors solidity, transition and melt with a total of 13 parameters by the following equations:



The coefficients  $b_5$  and  $b_6$  mark the seal- or freezing line at which material freezes. The model was originally developed for saline solutions. However, in the present expansion, it describes the measurement data of thermoplastics, elastomers and thermosets quite good as well and is used in well-known injection molding simulation software like Mouldflow, Moldex, CadMould and Sigmasoft.

$$\mathbf{Vsp(P,Tt) = Vsp(0,T) [ 1 - C \ln( 1 + P / B(T) ) ] + Vt(P,T) ; Tt(P) = b_5 + b_6 (P)}$$


$T < Tt(P)$  (Solid area):

$$\begin{aligned} Vsp(0,T) &= b_{1s} + b_{2s} (T - b_5) \\ B(T) &= b_{3s} \exp( - b_{4s} (T - b_5) ) \\ Vt(P,T) &= b_7 \exp( b_8 (T - b_5) - b_9 (P) ) \end{aligned}$$

$T > Tt(P)$  (Melting area):

$$\begin{aligned} Vsp(0,T) &= b_{1m} + b_{2m} (T - b_5) \\ B(T) &= b_{3m} \exp( - b_{4m} (T - b_5) ) \\ Vt(P,T) &= 0 \end{aligned}$$

# THIS IS RHEOLOGY



**GOETTFERT**  
THIS IS RHEOLOGY

**GOETTFERT Inc.**

Rock Hill, SC 29730  
USA

☎ +1 803 324 3883

✉ info@goettfert.com

**GÖTTFERT**  
THIS IS RHEOLOGY

**GÖTTFERT | Werkstoff-  
Prüfmaschinen GmbH**

74722 Buchen

☎ +49 (0) 62 81 408-0

✉ info@goettfert.de

**GÖTTFERT**  
CHINA LIMITED

**GOETTFERT (China) Ltd.**

Beijing 100027  
CHINA

☎ +86 10 848 320 51

✉ info@goettfert-china.com