



service



GÖTTFERT



THIS IS RHEOLOGY

# MAINTENANCE AND CALIBRATION

## SERVICE BASIC\*

Maintenance and calibration of the instrument with traceable measuring devices.

*After successful completion a maintenance check list (inspection protocol) will be provided.*

\* SERVICE BASIC is no longer provided for manually operated melt flow indexers of the current model series.

## SERVICE BASIC PLUS

Maintenance with factory calibration protocol: Maintenance and calibration of the instrument with traceable measuring devices, including statistical measurement uncertainty, which has been determined by an evaluation of the statistical measurands on several different instruments.

*After successful completion a factory calibration protocol (inspection protocol) will be provided.*

## SERVICE ISO 17025

Maintenance and calibration of the instrument with traceable measuring devices, including an indication of the lowest measurement uncertainty, which has been determined directly at the test device on site.

*After successful completion a verification report will be provided according ISO 17025.*

*Further information can be found on our homepage:*

[www.goettfert.com/service/maintenance](http://www.goettfert.com/service/maintenance)

